

PACBIO® GUIDELINES FOR SUCCESSFUL SMRTbell™ LIBRARIES

I. Sample requirements for PacBio sequencing

The Pacific Biosciences® library preparation process does not utilize amplification techniques and resulting library molecules are directly used as templates for the sequencing process. As such, the quality of the DNA starting material will be directly reflected in the sequencing results. Any irreversible DNA damage present in the input material (e.g., interstrand crosslinks, etc.) will result in impaired performance in the system. High-quality, high-molecular-weight genomic DNA is imperative for obtaining long read lengths and optimal sequencing performance.

II. General guidelines for handling high-molecular-weight DNA

In general, the following precautions need to be taken when handling DNA:

- Avoid overdrying of genomic DNA. Allow the DNA to air dry. Do not heat when drying in a speed-vac.
- DNA should be eluted in neutral, buffered solution (e.g., 10 mM Tris Acetate or Tris-HCl, pH 8) and stored in TE (10 mM Tris, pH 8, 1mM EDTA)*. Avoid eluting in RNase-free H₂O or unbuffered solutions.
- PCR products should be clean amplicons, without non-specific products or multiple bands.
- If gel purification is required, avoid using ethidium/UV based visualization methods. One alternative is to use SYBR® Safe (Invitrogen) and visualize with blue light.
- To help resuspend the DNA, carefully invert the tube several times after adding buffer and/or tap the tube gently.
- Alternatively, allow the DNA to stand in buffer overnight at 25°C to resuspend.
- Overheating can introduce DNA damage. Inactivate DNAase as recommended by the vendor kit. It is best to avoid heat inactivation when possible. An alternative is AMPure® purification.
- Avoid vortexing genomic DNA when possible as vortexing can cause shearing of the DNA.
- DNA storage conditions: 4°C (short-term); –20°C / –80°C (long-term).
- Repeated freezing and thawing of genomic DNA should be avoided as this will lead to DNA shearing.

*Note: EDTA must be removed prior to library preparation. This can be achieved during the initial AMPure purification.

III. Important measures impacting DNA quality

To maximize read length and quality, it is **essential** that the DNA sample:

- is double-stranded; single-stranded DNA cannot be used to generate the sequencing template.
- has not undergone multiple freeze-thaw cycles as they can lead to DNA damage.
- has not been exposed to high temperatures (e.g. > 65° C for 1 hour) or pH extremes (< 6 or > 9).
- has an OD₂₆₀/OD₂₈₀ ratio of 1.8 to 2.0.
- has an OD₂₆₀/OD₂₃₀ ratio of ~2.0.
- does not contain insoluble material.
- does not contain RNA contamination.
- has not been exposed to intercalating fluorescent dyes or ultraviolet radiation. SYBR dyes are not DNA damaging, but avoid ethidium bromide.
- does not contain denaturants (e.g., guanidinium salts or phenol) or detergents (e.g., SDS or Triton-X100).
- does not contain carryover contamination from the original organism/tissue (e.g., heme, humic acid, polyphenols, etc.)

IV. DNA sample quality assessment:

A thorough DNA quality check is required prior to submitting DNA for PacBio sequencing. The following recommendations to ascertain DNA integrity, purity, and concentration are recommended:

- 1) **Gel images of DNA sample:** Genomic DNA integrity can be assessed by agarose gel electrophoresis. For best results, DNA samples must show no signs of degradation, which is evidenced by smeared DNA bands. The presence of one predominant band showing high MW DNA with no degradation is optimal. A good practice is to indicate relevant marker sizes, and the amount of sample loaded in the agarose gel. For amplicon or cDNA samples, a Bioanalyzer® trace can be used as an alternative.

If gel purification of your DNA sample is required, we recommend using SYBR Gold or SYBR Safe coupled with blue light for visualization. Do not use ethidium bromide and/or UV light since they can induce DNA damage. The SYBR stains can be easily removed from nucleic acids during the gel-extraction process, e.g., using Qiagen® gel-extraction kits.

- 2) **Purity of your DNA sample:** DNA purity can be determined by using the NanoDrop® instrument or other spectrophotometers. Readings of both A260:A280 and A260:A230 ratios need to be obtained:

260/280:

The ratio of absorbance at 260 nm and 280 nm is used to assess the purity of DNA. A ratio of ~1.8 is generally accepted as “pure” for DNA, but is dependent on the nucleotide composition of the submitted sample.

A low A260/A280 ratio may indicate the presence of protein, phenol, or other contaminants that absorb strongly at or near 280 nm. Sometimes it may be caused by a very low concentration of nucleic acid.

High 260/280 ratios are not indicative of an issue.

Ensure DNA measurements are conducted in a buffered environment such as (TE or Tris HCl, pH8). Measurements are sensitive to small changes in the pH of the solution which will cause the 260/280 ratio to vary. Acidic solutions will skew the 260/280 ratio lower, while basic solutions will skew the ratio higher.

260/230:

The 260/230 ratios provides a secondary measurement of DNA purity to make inferences about the quality of sample extraction. Readings to determine purity are often higher than the respective 260/280 values. Expected 260/230 values are commonly in the range of 2.0-2.2. Abnormal 260/230 values may indicate a problem with the sample extraction procedure.

A low A260/A230 ratio may be the result of:

- Carbohydrate carryover (often a problem with plants).
- Residual phenol from nucleic acid extraction.
- Residual guanidine (often used in column-based kits).
- Glycogen used for precipitation.

A high A260/A230 ratio may be the result of:

- Making a blank measurement on a dirty pedestal of a Nanodrop instrument.
- Using an inappropriate solution for the blank measurement.

The blank solution should be the same pH and of a similar ionic strength as the sample solution.

- 3) **Concentration of your DNA sample:** Accurate quantitation of DNA concentration is critical for the PacBio® template preparation procedures. Traditional spectrophotometric assays cannot determine DNA concentrations <5 ng/μl. More importantly, almost all spectrophotometric assays do not distinguish between different types of

nucleotides (e.g., double-stranded DNA, RNA, dNTPs, and single-stranded DNA). Therefore, while the presence of single-stranded DNA will not impair library preparation, this will result in inaccurate yield quantitation.

For PacBio library preparation, it is critical to determine the concentration of the double-stranded DNA, since only double-stranded DNA will be converted into sequencing templates. RNA, dNTPs, and single-stranded DNA included in the concentration measurement will skew the concentration reading. Therefore, it is highly recommended to use the PicoGreen® assay or a Qubit® fluorimeter for quantitation purposes.

Please closely follow the recommended guidelines provided by the respective vendors when carrying out quantitation assays. In particular, the following steps should be observed: 1) periodically have the instrument calibrated (preferably by the vendor); 2) conduct a standard curve alongside samples when using the PicoGreen or Qubit assays; 3) perform replicate readings of concentration and use the average of replicates as the final concentration whenever possible. Accurate, consistent pipetting skill is needed to obtain reliable quantitation information.

V. Estimating Library yield on the PacBio® RS II (Table 1)

The PacBio Sample Calculator should be used to estimate yield for all samples. The table below provides estimates of expected yields for various DNA libraries. Please note the assumptions used to generate the table.

Library Insert Size*	Recommended Quantity for submission	Min Concentration Required (Post-Shearing)	Est. Total Yield (Range)	
			MIN	MAX
250 bp	600 ng	250 ng	60 GB	125 GB
500 bp	600 ng	250 ng	10 GB	20 GB
1 kb	1.2 µg	500 ng	90 GB	180 GB
2 kb	1.2 µg	500 ng	45 GB	90 GB
5 kb	2.4 µg	1 µg	45 GB	91 GB
10 kb	2.4 µg	1 µg	20 GB	45 GB
10 kb (AMPure kit)	10 µg	5 µg	90 GB	182 GB
20 kb (AMPure kit)	15 µg	5 µg	45 GB	91 GB
20 kb (BluePippin™ kit)	15 µg	5 µg	9 GB	18 GB

* Amounts recommended for submission represent quantities needed for one SMRTbell library prep and includes extra quantity needed for any additional QC and conservative excess. Reported library yield is based on an assumption of a DNA loading concentration of 50 ng/µl and throughput of 200 MB per SMRT Cell using P4-C2 chemistry. For insert sizes ≥1 KB, a magnetic bead loading protocol is used in the SMRT Cell calculation. Two size-selection protocols for large-insert libraries are available using either AMPure or BluePippin strategies. Actual results may vary.

VI. Accepted Volume and Buffers

- The maximum volume should not exceed 130 µl as this is the upper limit for some of the DNA shearing protocol. The upper limit for 20kb G-Tube® shearing is 200 µl.
- When resuspending DNA prior to generating SMRTbell libraries, such as during AMPure purification, avoid buffers containing EDTA to prevent enzymatic inhibition during downstream sample library preparation.
- Upon receipt of sample, performing AMPure purification to transfer DNA into a stable and appropriate buffer for downstream processing of PacBio SMRTbell templates is recommended. The additional benefit of this step is to further purify genomic DNA and remove carry-over contaminants.
- DNA can be dissolved in Tris buffer (e.g., 10 mM Tris, pH 7.0 – pH 8.0). Do not use nuclease-free water as this is insufficient for long-term DNA stabilization.

VII. Guidelines for SMRT® library preparation with high molecular weight and clean DNA

These are general recommendations to help obtain high molecular weight DNA.

Listed third party products are not officially endorsed by PacBio and are only provided as possible options.

1. Before DNA extraction:
 - a. Avoid incubation in complex or rich media.
 - b. Harvesting from several replicate cultures rather than a single, high-density culture during early- to mid-logarithmic growth phase is preferred.
 - c. Extraction of small volumes is preferred over large volumes to avoid accumulating high concentrations of potentially inhibiting secondary components.
2. Options for DNA Extraction:
 - a. Qiagen® MagAttract® HMW kit (100-200 kb) ([Product Information](#))
 - b. Qiagen Genomic-tip kit (50-100 kb) ([Product Information](#))
 - c. Qiagen Gentra® Puregene® kit (100-200 kb) ([Product Information](#))
 - d. Phenol-chloroform extraction (PacBio [SampleNet Protocol](#))
 - i. Ensure phenol is fresh and not oxidized; use within three months of opening of reagent bottle.
3. DNA cleanup before library prep
 - a. Purification of DNA with AMPure kit (Default recommendation, for all users)
 - b. Mo Bio® kit (MOBIO PowerClean® 50 prep kit, for highly contaminated samples)([Product Information](#))
 - i. Protocol requires approximately 30-45 min. Multiple samples can be prepped simultaneously.
 - ii. Some modifications to the protocol are suggested to maintain high-molecular-weight DNA and minimize damage:
 1. After adding Buffer 2, quickly vortex and add Buffer 3. Minimize time exposure of DNA in Buffer 2 to prevent damage.
 2. After adding Buffer 3 and vortexing, add 1 µl of glycogen and vortex before incubation on ice.
 3. During the final elution step, elute at 50 µl, spin for a second in the microcentrifuge, and incubate for 1 min at room temp. Proceed as stated in the protocol by adding 50 µl of elutant to the membrane for a minute incubation and final spin of 2 minutes.

WARNING: DNA recovery is low following this procedure; use as necessary. A 10 µg genomic DNA sample per column resulted in 30-50% recovery. A lower sample input will result in higher recovery, whereas a high sample input results in lower recovery.

4. Shearing:
 - a. G-Tube purification is recommended and is a preferred method due to the ease of use.
 - b. Hydroshear is alternative which allows shearing up to 20 kb. Please note use of Hydroshear requires maintenance to prevent frequent clogging of samples.
 - c. Covaris® E220 for fragments <5 kb.

For Research Use Only. Not for use in diagnostic procedures. © Copyright 2013, Pacific Biosciences of California, Inc. All rights reserved. Information in this document is subject to change without notice. Pacific Biosciences assumes no responsibility for any errors or omissions in this document. Certain notices, terms, conditions and/or use restrictions may pertain to your use of Pacific Biosciences products and/or third party products. Please refer to the applicable Pacific Biosciences Terms and Conditions of Sale and to the applicable license terms at <http://www.pacificbiosciences.com/licenses.html>.

Pacific Biosciences, the Pacific Biosciences logo, PacBio, SMRT and SMRTbell are trademarks of Pacific Biosciences in the United States and/or certain other countries. All other trademarks are the sole property of their respective owners.