

## **SOP for Thawing, Culturing, and Freezing Human iPSCs Stanford CVI Biobank**

The purpose of this document is to outline the standard operating procedures (SOP) for thawing and culturing human iPSCs.

### **From Life Technologies:**

Essential 8™ (E8) Medium: A1517001

### **From Stem Cell Tech:**

Gentle Cell Dissociation Reagent (EDTA): 07174

### **From Fisher Scientific:**

mTeSR: NC0722499

DMEM/F12 Medium: 10-092-CM

Y27632 2HCl (ROCK Inhibitor) 50mg: 50-863-7

BD Matrigel, hESC Qualified: 354277

Bambanker: NC9582225

**Note: This is for thawing one vial of iPSC**

### **Thaw Protocol:**

1. Dilute Matrigel 1:200 in DMEM/F12 and coat wells of a 6-well plate. Incubate the matrigel coated plates for 30 minutes at 37°C.
2. Prepare 2 ml of mTeSR medium in a 15 ml tube for washing.
3. Prepare another 2 ml of mTeSR medium + 10uM ROCK Inhibitor in a 15 ml tube for resuspension.
4. Remove frozen vial of iPSCs from liquid nitrogen storage and place in 37°C water bath until most, but not all, cells are thawed.
5. Remove cells from the freezing vial and add to the prepared tube of washing mTeSR medium.
6. Centrifuge at 200 g for 5 minutes and discard the washing supernatant.
7. Resuspend the cells in 2 ml of mTeSR medium plus ROCK inhibitor.
8. Aspirate excess Matrigel from well surfaces and plate cells into 1 well of a six well plate. Incubate cells in a 37°C, 5% CO2 incubator.
9. The next day, replace the spent mTeSR + 10uM ROCK Inhibitor medium with 1 ml each of fresh mTeSR and E8 medium (no ROCK Inhibitor in either).
10. The next day, replace with only E8 media. Replace the medium daily until the cells reach 85% confluence.

### **Passage of hiPSC with EDTA**

1. Ideally cells should have reached 85% confluence in 4 days. Use a split ratio of 1:6.
2. Aspirate culture medium

3. Add 1 mL per well of Gentle Cell Dissociation Reagent (EDTA), incubate for 7 min at RT (in hood)
4. During this time set up a 15 mL Falcon tube with fresh E8-RI medium (i.e. 12 mL for 1:6 split)
5. Aspirate EDTA from well
6. With a P1000 tip, add 1 mL of the E8-RI medium and blast against cell surface to dissociate cells. Cells should come off easily. Transfer to the Falcon tube.
7. Matrigel-coated plate: aspirate excess Matrigel
8. Plate out cells at 2 mL per well.
9. The next day, replace the spent E8 + ROCK Inhibitor medium with fresh E8 medium (no ROCK Inhibitor). Replace the medium daily until the cells reach 85% confluence.

### **iPSC Freezing Protocol**

1. Aspirate culture medium
2. Add 1 mL of EDTA per well, incubate for 7 min at RT (in hood)
3. Aspirate EDTA from well
4. With a P1000 tip, add 1 mL of freezing medium (Bambanker) and blast against cell surface to dissociate cells. Cells should come off easily. Transfer cell suspension into freezing tube (1 well per tube).
5. Place freezing tubes in a freezing container and place the container in -80°C freezer overnight.
6. After 24 hours move the freezing vials to liquid nitrogen storage.

### **For more information please contact:**

#### **Joseph C. Wu, MD, PhD**

Director, Stanford Cardiovascular Institute  
Professor, Dept of Medicine & Radiology  
Email: [joewu@stanford.edu](mailto:joewu@stanford.edu)

#### **Ioannis Karakikes, PhD**

Co-Director, SCVI Biobank  
Email: [ioannis1@stanford.edu](mailto:ioannis1@stanford.edu)

#### **Justin Vincent**

Manager, SCVI Biobank  
Email: [justin81@stanford.edu](mailto:justin81@stanford.edu)

#### **Rinkal Chaudhary**

Co-manager, SCVI Biobank  
Email: [rinkalc@stanford.edu](mailto:rinkalc@stanford.edu)