



Automated Intradural Tumor Segmentation and Sub-Type Classification in Spinal MRI Scans

Golnaz Moallem, Ph.D.

Postdoc Research Fellow, PIMed Lab

Department of Radiology, Stanford School of Medicine

January 2024

Outlines

➤ Introduction

- What are intradural spinal tumors?

➤ Data

➤ Methods

- Baseline framework
- Proposed framework


➤ Results

➤ Challenges and Next Steps


Introduction: What are Intradural Spinal Tumors?

Spinal Tumors

- Large heterogenous group of abnormal mass of tissue affecting the spinal cord
- Developed within or near the spinal cord
 - **Intradural:** growth within the dura of the spinal cord
 - Extradural: mass within the spinal column
- Malignant or Benign
- Ramifications: pain, neurological damage, loss of mobility, etc.




Intradural Spinal Tumor Treatment Workflow



Intradural Spinal Tumor Detection Framework

Tumor Detection Supplemental Tool



Data

Data

- **84 intradural spinal tumor cases** from the Stanford University hospital database
 - Surgical cases; pathology images available
 - T1-weighted, T2-weighted, and post-contrast T1-weighted MRI scans
 - The pathology slides are used to determine the tumor sub-types


Sagittal View



Axial View




Pathology Slide



Data

- Variation in the field of view; different spinal section
 - Cervical: 29 cases
 - Thoracic: 31 cases
 - Lumbar: 24 cases

Cervical case



Post-Contrast T1 scan

Thoracic case



Post-Contrast T1 scan

Lumbar case




Post-Contrast T1 scan

Data

- Variation in the appearance of the tumor; different tumor sub-types
 - Ependymoma: 21 cases
 - Meningioma: 23 cases
 - Schwannoma: 26 cases
 - Others (Hemangioblastoma, Neurofibroma, etc.): 16 cases

Ependymoma




Post-Contrast T1 scan

Meningioma




Post-Contrast T1 scan

Schwannoma



Post-Contrast T1 scan

Others



Post-Contrast T1 scan


Data

- Variation in the volume intensity distributions; different manufacturers
 - GE MEDICAL SYSTEMS: 68 cases
 - SIEMENS: 10 cases
 - Others (TOSHIBA, PHILIPS , HITACHI): 6 cases

GE MEDICAL SYSTEMS



SIEMENS



HITACHI



Methods

Methods: Baseline Framework

Prediction Annotation


Post-Contrast T1-Weighted Spinal
MRI Scan



3D U-Net Tumor Segmentation
Model



Segmented Tumor



Methods: Proposed Framework

Prediction Annotation

T1-weighted scan



Step 1:
Spinal Canal
Segmentation

Spinal Canal Mask




Post-contrast T1-weighted scan



Step 2:
Intradural Tumor
Segmentation

Segmented Tumor




Methods: Proposed Framework - Spinal Canal Segmentation

T1-weighted scan



Spinal Canal Segmentation
Attention U-Net Model



Spinal Canal Mask



Methods: Proposed Framework –Tumor Segmentation


Post-contrast T1-weighted scan



Spinal Canal Mask




Intradural Tumor Segmentation



Segmented Tumor



Swin UNETR Segmentation Model [1]




Prediction Annotation

[1]: Hatamizadeh, Ali, et al. "Swin unetr: Swin transformers for semantic segmentation of brain tumors in mri images." International MICCAI Brainlesion Workshop. Cham: Springer International Publishing, 2021.

Methods: Proposed Framework - Data Annotation Process

- **Tumor Segmentation:** Manually annotated the tumors using 3D Slicer tool[1] which were later verified by the spinal neurosurgical expert
- **Spinal Canal Segmentation:** Developed a unique algorithm which automatically segmented the canal for a group of training examples




Results


Segmentation Example

 Prediction  Annotation

Post-contrast T1-weighted scan




Masked post-contrast T1-weighted scan




Segmentation Example

 Prediction  Annotation

Post-contrast T1-weighted scan



Masked post-contrast T1-weighted scan




DICE = 0.78


Segmentation Example

 Prediction  Annotation

Post-contrast T1-weighted scan



Masked post-contrast T1-weighted scan



DICE = 0.65

Quantitative Results: Spinal Canal Segmentation

Spinal Canal Segmentation			
	No. of testing cases	Patient-level Sensitivity	Mean DICE
Attention UNET-based model	150 (normal, tumor, fusion, cyst, osteotomy)	95%	0.85

Quantitative Results: Intradural Tumor Segmentation

Intradural Tumor Segmentation				
K-Fold Cross Validation (K = 4), No. of Tumor Cases: 84 cases				
Fold	Patient-level Sensitivity		Mean DICE	
	Baseline	Proposed	Baseline	Proposed
Fold 1	67%	80.9%	0.47	0.75
Fold 2	75%	90.5%	0.52	0.65
Fold 3	63%	85.7%	0.41	0.63
Fold 4	62%	82.6%	0.53	0.62
Average	66.8%	84.9%	0.48	0.66

Intradural Tumor Segmentation	
No. of Non-Tumor Cases (31 cases)	
Patient-level Specificity	
Baseline	Proposed
86%	92%

Challenges and Next Steps

Challenges & Next Steps

- **Limited tumor sub-type categories**

Employ semi-supervised techniques for tumor sub-type classification

- **Multiple manufacturers and variation in intensity distributions**

Develop an intensity histogram alignment to standardize the intensity distributions

Conclusion

- An AI-powered intradural spinal tumor diagnosis tool presented to assist the radiologist/surgeon in identifying these tumors in MRI scans
- For this purpose, an automated segmentation model proposed for segmenting intradural spinal tumors in MRI scans
- Quantitative and qualitative results for this segmentation model demonstrated on a cohort of surgical intradural tumor cases

Thank You!



Mirabela Rusu, PhD

Assistant Professor Of Radiology and, By
Courtesy, Of Urology



Atman Desai, MD

Spine Neurosurgeon
Clinical Associate Professor, Neurosurgery



Greg Zaharchuk, MD

Professor of Radiology